

## Sher™ Comfort Series

The Sher™ comfort Series offers a full range of half-sole inserts to address a variety of conditions. All made from high-quality EVA/fabric laminate with different firmnesses to meet your patients' needs. Dual-layer construction for durable, long term-wear.



### Model 717 Softest 35-40 Durometer

- Soft, gentle cushioning with resilient support
- Full contact with arch and heel for total comfort
- General use without modifications, trimmable

Ideal for: geriatric foot support, general comfort, plantar fasciitis



### Model 718 Soft 35-40 Durometer with Poron Heel Pad

- Gentle cushioning support for total comfort
- Energy-absorbing Poron heel pad softens heel strike
- General use without modifications, trimmable

Ideal for: plantar fasciitis, heel pain, sports and activities



### Model 719 Medium/Firm 55-60 Durometer

- Firmer density support against the foot and arch
- Full contact with arch and heel for total comfort
- General use without modifications, trimmable

Ideal for: plantar fasciitis, arch pain and support



### Model 720 Dual Density Top 35-40, Bottom 20-25 Durometer

- Dual density construction with firmer support against the foot
- Softer bottom for cushioning support
- Conforms to the foot with wear
- Provides comfort to high pressure areas
- Superior shock absorption and cushioning
- General use without modifications, trimmable

Ideal for: plantar fasciitis, foot pain, sports, arch pain and support



### Model 721 Firmest 75-80 Durometer

- Compression resistant, retains shape and form after extended wear
- Sturdy and long lasting
- Can be modified, grindable for a customized shape
- Molded heel wedge
- Excellent base for custom orthotics

Ideal for: excessive pronation, arch support, foot stabilization

## Deluxe Nylon "E-Z" Post™

491 - 3/4 Length  
493 - Full Length

The combination of one of our best selling "Off the Shelf" inlays with our newly introduced "E-Z" Post™.

Available in both 3/4 and full length.

A 1/16" gray Poron center laminated to blue ProCell.

Laminates are attached to our semi-rigid polypro shell.

The "E-Z" Post semi-rigid shell has a simple snap-in ability to accurately post the shell at 2° (blue), 4° (red) or Neutral (gray). The plates reverse right to left for either a Medial or Lateral Post.

The orthotic may also be custom molded using a heat gun.

Available in Women's 5 - 12 1/2

Men's 7 - 13

Also available in leather 490 3/4 Length and the 492 Full Length



Call Your Local Representative:

**SROUFE**  
HEALTHCARE PRODUCTS  
601 Sroufe Street • P.O. Box 347 • Ligonier, IN 46767  
TEL (888) 894-4171 • FAX (877) 894-4092  
www.sroufe.com



Professional Medical  
Grade Orthotics

Sher Brand Orthotics  
are made from the finest  
quality materials in time  
tested designs based on  
hand craftsmanship.

Discover the Quality Sroufe  
Healthcare Products lines.



## Customize

Right from your Business Card, Letterhead or Graphic Logo Files.



**460 - 3/4 Length Leather DS**  
**462 - Full Length Leather DS**

### Premium Leather Top Cover

Supple, breathable leather conforms to the foot with use for even greater comfort. Helps prevent and relieve calluses on the forefoot and toes.

### High Density Poron Cushion Foam

Durable foam backing has high rebound and is resilient to long term wear.

### Semi-Rigid Plastic Arch Support

**Shells** - Gives arch support to help relieve pain and stress. Semi-rigid polypropylene shells stabilize the foot, controlling excessive motion and rotation to reduce the stress on the knee and hip joints. Shells can be heat molded using a heat gun.

- Promote your business with a top quality product and retain customers for repeat orders.



**481 - 3/4 Length Deluxe Nylon**



**483 - Full Length Deluxe Nylon**

**Nylon Top Cover** - Breathable nylon minimizes friction, shear and heat build up to help prevent and relieve blisters and calluses.

**High-Density Poron and ProCell Cushion Foams** - Durable foam backing has high rebound and is resilient to long term wear.

- Gives arch support to help relieve pain and stress. Semi-rigid polypropylene shells stabilize the foot, controlling excessive motion and rotation to reduce the stress on the knee and hip joints. Non skid tear drop at the heel holds orthotic in place. Shells can be heat customized and molded using a heat gun.



**4200 Tri-Lam**



**4210 Bi-Lam**

**Tri-Laminate Construction** - Provides greater relief and comfort to diabetic feet with cushioning Poron middle layer.

**Heat Moldable** - Provides a total contact, customized fit to relieve pressure points.

**Deep Heel Cup** - for increased heel and foot stabilization

**Moderate longitudinal arch support**

**Easy to trim and fit in most diabetic shoes**



**4100 - PF Night Splint**

**Adjustable Flexion Straps** - Controls the amount of dorsi-flexion.

**Padded Straps with Lock Buckles** - Comfortable wear and easy to apply.

**Lightweight Shell and Padded Liner**  
**Non-Skid Pad**



**4102 - Anterior Night Splint**

**Lightweight and Cool Low Profile Design** - Easier for patients to wear at night. less bulk and less heat build-up.

**Padded Non-Stretch Straps** - Holds the splint in place to keep proper foot position to prevent passive plantar flexion.



**1765 SAS Sher Ankle Stabilizing Brace**

**Stabilizing Figure-Eight Straps** - Mimic ankle taping, locks heel and foot in place.

**Wide Elastic Cuff Closure** - Holds straps in place and provides added stability and compression.

**Lace Up Design** - Easy to apply and gives a conforming fit

**Ballistic Nylon Material** - Extremely lightweight and thin with long-lasting strength.

**Low Profile** - Easily to wear and fits in any type of shoe.

**Breathable Tongue and Heel** - Advanced material wicks moisture away from the skin and allows air to circulate. Tongue is sewn in place to reduce brace movement.

PDAC Approved for HCPCS Code L1902



**501 "Sher Slim Line"**

Designed with women's feet in mind. A slimmer profile to fit into women's dress shoes. Made using our semi-rigid polypro shell covered with a grey poron laminated to leatherette.

Sher



601 Sroufe Street  
P.O. Box 347  
Ligonier, IN 46767

Phone: (888)-894-4171  
Fax: (260)-894-4092  
[www.shercoinc.com](http://www.shercoinc.com)

### Sadmerc Coded:

Our diabetic inlays are Sadmerc coded for reimbursement using 5512 that includes the 5512, 4200 and the 4210. The 5513B is coded for reimbursement using 5513.

### Custom Personalization:

Call us about personalizing your order with your logo. Minimum of 48 pairs - mixed sizes. Available in 2 - 3 weeks on the initial order and 7 - 10 days on future orders.

Visit our website at [www.shercoinc.com](http://www.shercoinc.com)